

Life Support Systems offers multiple options for oxygen systems for patients who need supplemental oxygen while at home or on the go. We respond to the needs of every patient individually by providing the oxygen system most appropriate for his or her lifestyle. Oxygen concentrators, stationary and portable liquid oxygen systems, and high pressure cylinders are provided according to the physician's prescription and with consideration for the activity level of each patient.

Oxygen concentrators allow patients to receive uninterrupted continuous oxygen without running empty. Concentrators are more mobile than stationary liquid reservoirs or oxygen cylinders, and can be wheeled to any location within the home and within a reasonable distance outside the home.

Oxygen systems (portable or stationary) may be appropriate for patients who want to remain active and participate in exercise rehabilitation, or for patients who have high continuous liter flow requirements.

Cylinder oxygen systems are for patients who require a low liter flow or high liter flow for intermittent usage.

### Oxygen Equipment

Airsep Lifestyle Power Pak	<a href="#">Manual</a>
Caire Lib 10 20 45 60	<a href="#">Manual</a>
Companion 1000	<a href="#">Manual</a>
Easy Mate	<a href="#">Manual</a>
Eclipse Travel Concentrator	<a href="#">Manual</a>
Everflow Concentrator	<a href="#">Manual</a>
Helios	<a href="#">Manual</a>
Millennium Concentrator	<a href="#">Manual</a>
Perfecto 2 Series Concentrator	<a href="#">Manual</a>
Platinum Series 5 & 10 Concentrator	<a href="#">Manual</a>
Spirit 300	<a href="#">Manual</a>

